



TESTQUAL

PROFICIENCY TESTING SCHEMES

TestQual, S.L.

(Proficiency Testing Schemes)

*Pol. Industrial Oeste
Av. Principal, Parcela 21/1
CP 30169, San Ginés, Murcia
Telephone: 868 949 486 / 676 367 555*



TestQual 146 A+ B PROTOCOL

Pesticides Residues in Grapefruit and Peach

0. GLOSARY AND ABBREVIATIONS

Text	Abbreviation
TestQual	TQ
Proficiency test	PT / P.T.
Limit Of Quantification	LOQ
NA	Not Analysed

1. INTRODUCTION

This document describes the **protocol** of the **TestQual 146 A+ B** Proficiency Test (P.T.), belonging to the analysis of **pesticides** in **Grapefruit and Peach**.

In the present document is detailed how to start working with TestQual, send your application to participate, the statistic that will be applied and information about the evaluation report.

The following proficiency scheme is planned, scheduled and will proceed as one proficiency test for two different matrixes that for technical reasons on our website are separated in two: TestQual 146A Pesticides in grapefruit and TestQual 146B pesticide in Peach.

There is a special price for requesting both matrixes, you can ask for your quote using the contact data in section 11 of the present protocol.

TestQual, S.L. is committed to maintaining confidentiality with the information of each laboratory from the beginning of the proficiency test.

2. OBJECTIVE

The objective of the **TestQual 146 A+ B** Proficiency Test is to evaluate the quality and accuracy of the results sent by the participating laboratories. Because of this, proficiency testing is an essential element of laboratory quality assurance. It will help to control and detect errors in their results or methods of analysis.

3. CALENDAR

The following table shows the program for this proficiency scheme (for both samples):

Date	Activity	Carried out by
21/May/2021 (Week 20)	Final date to submit applications	Participants
15/Jun/2021 (Week 24)	Sample delivery	TestQual
09/Jul/2021 (Week 27)	Final date to submit results	Participants
30/Jul/2021 (Week 30)	Final report (Email and/or client area)	TestQual

The dates of this calendar might be slightly changed according to the development of the proficiency test during the year. However, any change would be notified to all participants announcing it on our website www.TestQual.com.

The **coordinator** of this proficiency test will be Jose Pedro Navarro. Any question regarding the development of the proficiency test can be consulted by email to jpnavarro@testqual.com.

4. REGISTER AND PARTICIPATION REQUEST (APPLICATION FORM)

NEW CLIENT

If your laboratory has not participated before in one of our proficiency tests you will have to register on the [REGISTER](#) form.

Once you have completed and sent the form you will have to wait until the activation of the account from the website administrator. If some more information is needed someone from our team will get in contact with you through the phone or email you used during your registration.

You can find our contact data at the end of this protocol.

For those laboratories that require more than one contact per account or are in a situation not contemplated in this protocol will have to contact the organizer using the Contact tab to be instructed how to proceed.

APPLICATION FOR THE PROFICIENCY TEST

To participate in this proficiency test or be updated via email of any changes regarding this proficiency test is required to have a laboratory code. To get your laboratory code is needed to apply through the website or have it arranged by the coordinator at least 24h before the shipping of the samples.

In the Proficiency Tests Tab on our website you will have to select the proficiency test you want to participate, by clicking its name or the shopping cart you will enter the page with general information and a summary of that proficiency test, there you can find the present document (the protocol) and the button to start the application.

If you did not log in before you will be requested to do so and then the website will require you to submit your Limit Of Quantification (LOQ) for the parameters you will study. The compounds left as NA (NOT ANALYSED) will NOT appear in the Results form and therefore will not be able to send results for that parameters through the website.

Once the application has been sent, as soon as possible, it will be checked by the website administrator and you will be sent an email with the participation code. This code will be just known only by the organizer and the participant and will be kept confidential at all times, even after the proficiency test finishes.

To participate in the proficiency scheme (TQ146A + TQ146B) you will have to send the application for both proficiency tests on our website.

You can check on the dashboard of your client area if an application you sent has been accepted or is still pending.

Just one application per exercise can be sent by each laboratory, being not allowed for a laboratory to participate with two different codes.

The applications of the laboratories will be studied and accepted in base of the quantification limits of the analytes of the proficiency test and if the logistics allow the sample shipping without risk of deterioration.

According to the experience, TestQual can anticipate that the number of participants of this proficiency test will be around 15-20, being 11 the minimum participants for the proficiency test to take place.

5. PREPARATION AND QUALITY CONTROLS

TestQual 146A+ TQ146B scheme is a proficiency test based in the analysis of **pesticides** in **Grapefruit and Peach** that has been spiked with pesticide **standards**. The material will be bought in an ecological shop in Murcia and analysed by a subcontracted laboratory that holds the standard UNE-EN ISO/IEC 17025 into force.

The material is cut in very small pieces, dropped into liquid N₂, once it is fully frozen, it is grounded into a fine powder, puree or juice, which will be spiked with a solution with the analytes of the PROFICIENCY TEST and poured into a homogenizer with controlled temperature to ensure complete homogeneity.

Once the lot of samples is ready they will be stored in a temperature-controlled freezer below -20°C until the dispatch of the samples.

Before the samples are distributed, for the assessment of the homogeneity of the lot of samples prepared, ten samples from the lot will be selected randomly and analysed in duplicate by TestQual's collaborator laboratory under repeatability conditions. If the mean concentration obtained from this study is not within the planned range, the participants will be informed and a new distribution day might be scheduled if another spiking is deemed as necessary.

For stability assessment purpose, three samples are analysed, in duplicate, before, during and at the end (once all laboratories have sent the results) of the proficiency test.

6. TEST MATERIAL AND SAMPLE SHIPMENT

The shipment of the test materials (both matrixes) will take place on the date shown in the calendar, to the address provided by each laboratory in their client area. Specific delivery dates can change from the scheduled dates of the calendar, but all changes will be announced both in the website and by mail to the registered laboratories.

This PT will consist in a single round in which will be sent a sample of approximately **100-150 g** of test material. The samples will be sent by courier service (MRW, FedEx, DHL or TNT, depending on the destination). The material will be sent in insulated box that ensure the temperature conditions of the package during the whole shipment. The transit will be 1, 2 or 3 days to the destination country, depending on the location of the receiving laboratory. These boxes will be provided with either dry ice or cold packs to keep the temperature.

The shipping costs are not included in the price displayed on the website, which can only be seen if you are registered and logged in. To get an approximation of the shipping costs you can get your quotation by using the contact data at the end of this protocol.

A second test material for both matrixes can be requested date if necessary. If the package and/or the sample arrived damaged, defective or not valid the participating laboratory will have to notify of this to the coordinator before two working days to get another sample.

Before the shipment, TestQual will send the instructions for storage and analysis via email and confirm the distribution date. You can request a paper copy to be attached to the package and/or TestQual might decide to include it in addition to have it sent by email.

From TestQual we encourage our participants to read it carefully and follow its instructions, as it can help to conserve correctly the sample and increase the reproducibility of the analysis.

You can request a digital copy of this document by letting us know through any communication channel.

7. CONCENTRATION RANGES, SIGMA OBJECTIVE AND ANALYTES

The range of concentration for the target analytes of this proficiency test might be between **10** and **200 µg/kg** approximately.

The **sigma objective ($\hat{\sigma}$)** which works in this scheme will be the **25 % of the assigned value**. This value has been chosen according to the experience of similar proficiency tests organized by TestQual.

The **possible pesticides** in the Grapefruit and Peach are from the list below:

2-Phenylphenol	Boscalid	Chlorpyrifos Methyl	Demeton-S-methyl	Disulfoton
3,5-Dichloroaniline	Brodifacoum	Chlorthion	sulfone	Ditalimfos
3-Hydroxy-carbofuran	Bromacil	Chlorthiophos	Desmetryn	Diuron
4,4-	Bromocyclen	Cyanazine	Dialifos	Dodine
Dichlorobenzophenone	Bromophos-ethyl	Cyazofamid	Diazinon	Emamectin benzoate
Abamectin	Bromophos	Cyfluthrin	Dicaphon	B1a
Acephate	Bromopropylate	Cymoxanil	Dichlofenthion	Endosulfan-alpha
Acetamiprid	Bromuconazole	Cypermethrin	Dichlormid	Endosulfan-beta
Acetochlor	Bupirimate	Cyproconazole	Dichlobenil	Endosulfan-sulfate
Aclonifen	Buprofezin	Clethodim	Diclobutrazol	Endrin
Acrinathrin	Butafenacil	Clofentezine	Dichlofluanid	EPN
Alachlor	Butamifos	Clomazone	Diclofop-methyl	Epoxiconazole
Aldicarb	Butoxycarboxim	Cloquintocet-mexyl	Dicloran	Etaconazole
Aldicarb sulfone	Butralin	Chlorfenson	Dicrotophos	Ethion
Aldicarb sulfoxide	Buturon	Chlorotoluron	Dieldrin	Ethoprophos
Aldrin	Cadusafos	Chloroxuron	Diethofencarb	Etoxazole
Antraquinone	Captan	Chlorpropham	Difenoconazole	Ethiofencarb
Atrazine	Carbaryl	Chlorsulfuron	Difenoxuron	Ethiofencarb -sulfone
Azaconazole	Carbendazim	Chlorthal-dimethyl	Diflubenzuron	Ethiofencarb -sulfoxide
Azinphos-ethyl	Carbophenothion	Clothianidin	Diflufenican	Etofenprox
Azinphos-methyl	Carbofuran	Coumaphos	Dimethenamid	Ethofumesate
Azoxystrobin	Chlorantraniliprole	Kresoxim-methyl	Dimethoate	Etrimfos
Benalaxyl	Chlorbromuron	Crimidine	Dimethomorph	Famoxadone
Bendiocarb	Chlorfenapyr	Cyanofenphos	Dimoxystrobin	Famphur (Famophos)
Benfluralin	Chlorfenvinphos	Cyanophos	Diniconazole	Fenarimol
Benfuresate	Chlormephos	Cycloxydim	Dioxacarb	Fenazaquin
Bentazone	Chloroneb	Cyprodinil	Dioxathion	Fenbuconazole
Bifenthrin	Chloropropylate	Deltamethrin	Diphenylamine	Fenbutatin oxide
Bitertanol	Chlorpyrifos	Demeton-S-methyl	Dipropetryn	Fenclorphos

TestQual 146 A+B Pesticides in grapefruit and peach. Protocol. Rev. 00
TestQual, S.L. Proficiency Testing Schemes

Fenhexamid	Hexaconazole	Metribuzin	Picoxystrobin	Spiroxamine
Fenitrothion	Hexaflumuron	Mevinphos	Piperonyl butoxide	Sulfotep
Fenoxycarb	Hexazinone	Myclobutanil	Pyraclostrobin	Sulprofos
Fenpropathrin	Hexythiazox	Molinate	Pyrazphos	Tebuconazole
Fenpropimorph	Imazalil	Monocrotophos	Pyridaben	Tebufenozide
Fenpyroximate	Imazamethabenz-	Monolinuron	Pyrifenox	Tebufenpyrad
Fensulfothion	methyl	Monuron	Pirimicarb	Tebupirimfos
Fenthion	Imidacloprid	Napropamide	Pirimicarb-desmethyl	Tecnazene
Phenthoate	Indoxacarb	Neburon	Pirimiphos-ethyl	Teflubenzuron
Fenuron	Iprobenfos	Nitenpyram	Pirimiphos-methyl	Tefluthrin
Fenvalerate	Iprodione	Nitrofen	Pyriproxyfen	Terbacil
Fipronil	Iprovalicarb	Nitrothal-isopropyl	pp-DDE	Terbufos
Flonicamid	Isazofos	Norflurazon	pp-TDE(DDD)	Terbumeton
Fluazifop-P-butyl	Isocarbophos	Nuarimol	Prochloraz	Terbutylazine
Fluchloralin	Isofenphos	Ofurace	Procymidone	Terbutryn
Flucythrinate	Isofenphos-methyl	Omethoate	Propham	Tetraconazole
Fludioxinil	Isoproturon	op-TDE (DDD)	Profenofos	Tetradifon
Flufenoxuron	Lambda-Cyhalothrin	Oxadiazon	Profluralin	Tetramethrin
Flumetralin	Lenacil	Oxadixyl	Promecarb	Tetrasul
Fluometuron	Leptophos	Oxamyl	Prometryn	Thiabendazole
Fluotrimazole	Linuron	Oxamyl-oxime	Propachlor	Thiacloprid
Fluquinconazole	Lufenuron	Oxydemeton-methyl	Propamocarb	Thiamethoxam
Flusilazole	Malaoxon	Oxyfluorfen	Propanil	Thiodicarb
Flutolanil	Malathion	Paclobutrazol	Propargite	Thiobencarb
Flutriafol	Mecarbam	Parathion	Propetamphos	Thiometon
Folpet	Mefenpyr-diethyl	Parathion-methyl	Propiconazole	Tolclofos-methyl
Fonofos	Mepanipyrim	Pebulate	Propyzamide	Triadimefon
Formothion	Mepronil	Penconazole	Propoxur	Triadimenol
Phosalone	Metalaxyl	Pendimethalin	Prosulfocarb	Triazophos
Phosphamidon	Metamitron	Pentachloroanisole	Prothiofos	Trichloronate
Phosmet	Metazachlor	Permethrin	Pyridafenthion	Tridemorph
Furalaxyl	Methacrifos		Pyrimethanil	Trifloxystrobin
Furathiocarb	Methamidophos	1,1-(2,2-	Quinalpho	Triflumuron
HCH-Alpha	Methidathion	dichloroethylidene)	Quinoxyfen	Trifluralin
HCH-Beta	Methomyl	bis(4-	Quintozene	Vinclozolin
HCH-Delta	Methoxychlor	methoxybenzene)	Rotenone	Yodofenfos
HCH-Gamma (lindane)	Methoxyfenozide		Simazine	Zoxamide
Heptachlor	Metobromuron	(methoxychlor	Simetryn	
Heptachlor-epoxide	Metolachlor	metabolite)	Spinosad A+D	
Heptenophos	Methoprotryne		Spirodiclofen	
Hexachlorobenzene	Metoxuron	Phenmedipham	Spiromesifen	

8. RESULTS EXPRESSION

Each participant laboratory must analyse the sample received according to their routine procedure, and fill up the RESULTS form of its client are of the website www.TestQual.com with just one value per analyte/parameter.

The results should be expressed in **µg/Kg**. The number of significant figures and the units are to be chosen by laboratories and will be displayed in the report as received through the website.

The method used for the analysis of each compound informed should be sent when filling up the results form.

The organizer should get the results before the fixed data of the scheme.

If you have any problem logging in to your client area or submitting your results you can contact the coordinator of the PT for guidance or help.

Once the results are sent you can check if they were correctly recorded by accessing the detailed information of this proficiency test, which can be accessed in your client area.

9. STATISTICAL EVALUATION

TestQual will develop the following statistical evaluation for each matrix of this Proficiency scheme:

TestQual considers as an **extreme outlier** any data which differs more than **50 %** of the average of all results reported by the laboratories, according to the Harmonize Protocol of the IUPAC. These extreme values will not be included in the calculation of the assigned value.

Once received all the results, TestQual evaluates the unimodality of all the values by Kernel test, being explained in the final report which is the followed procedure in case there is more than one distribution.

The **assigned value (X)** is determined using the robust average of the results considered valid for statistical computing (after eliminating the extreme outliers), according to the standard ISO 13528 into force.

The **standard uncertainty (u_x)** is calculated using robust statistics from the following formula:

$$u_x = s^*/\sqrt{p}$$

Being s^* the robust standard deviation of the data and p the number of results considered.

The following condition must be fulfilled in order to discard the contribution of the uncertainty:

$$u_x \leq 0,3 \hat{\sigma}$$

In case this condition is not fulfilled, the participants of the scheme will be informed in the report, and the uncertainty will have to be taking into account for the assigned value assessment.

The **standard deviation for proficiency assessment**, also named **target standard deviation, ($\hat{\sigma}$)**, comes from this formula:

$$\hat{\sigma} = b_i \cdot X$$

Being $b_i = \%_{DSRA} / 100$, and $\%_{DSRA}$ is the assigned relative standard deviation.

In this case, the assigned relative standard deviation is **25 %**. This value is fixed previously by the organizer based in the experience of TestQual organizing similar proficiency tests.

Proficiency assessment (z-score): This parameter shows the competence and accuracy of the laboratory. It is calculated using the following formula:

$$z = (x_i - X) / \hat{\sigma}$$

Where x_i is the value reported by the laboratories, X is the assigned value, and $\hat{\sigma}$ is the target standard deviation for each analyte.

The criteria for defining the z-score values are:

	$ z $	≤ 2	Satisfactory
$2 <$	$ z $	≤ 3	Questionable
	$ z $	> 3	Unsatisfactory

False negatives: Any analyte not reported in the results that is in the sample above the limit of quantification previously established to the proficiency test established by the organization (**10 $\mu\text{g/Kg}$**). TestQual assigns to all false negatives a result equal to half the laboratory limit of quantitation (LOQ/2).

False positives: Those analytes reported in the results, which is not present in the test material, and is reported by the participant at concentrations higher than the limit of quantification of the P.T. (**10 $\mu\text{g/Kg}$**).

Testing for sufficient homogeneity:

Once the samples are prepared ten of them will be chosen at random and sent to be analysed by TestQual's collaborator laboratory. Once received the results, a statistical evaluation will be performed, according to the IUPAC Harmonic Protocol.

The acceptance criterion to ensure that the randomly chosen samples are homogeneous is that the square of the estimated sampling standard deviation is below the critical value for accepting proper homogeneity:

$$S_{sam}^2 < c$$

In the first place to check the criterion, S_{sam}^2 which is the estimated sampling standard deviation, was calculated from:

$$S_{sam} = \left(\frac{V_s}{2} - S_{an} \right)$$

Firstly V_s is the variance of the sums S_i :

$$V_s = \sum \frac{(S_i - \bar{S})^2}{m - 1}$$

Where S_i was obtained from the addition of each duplicate result from the homogeneity; \bar{S} is the mean of all S_i and m is the number of samples (10 samples).

And secondly S_{an}^2 , which is the experimental estimate of analytical standard deviation, is obtained following the next formula:

$$S_{an}^2 = \frac{\sum D_i}{2m}$$

where D_i is the result of the subtraction of each pair of replicates from the homogeneity and m is the number of samples.

In second place to check the criterion for sufficient homogeneity the critical value c was obtained from:

$$c = F_1 \cdot \sigma_{all}^2 + F_2 \cdot S_{an}^2$$

Being F_1 and F_2 constants with values equal to 1.88 and 1.01 respectively for 10 samples. S_{an}^2 has already been calculated and σ_{all}^2 is obtained from:

$$\sigma_{all}^2 = (0.3 \cdot \hat{\sigma})^2$$

where $\hat{\sigma}$ is the target standard deviation, which is calculated with the formula:

$$\hat{\sigma} = 0.25 \cdot \bar{X}$$

Being \bar{X} , the mean of the 20 values from the homogeneity.

Testing for sufficient stability:

Three samples will be analysed, in duplicate, before, during and at the end (once all laboratories have sent the results) of the proficiency test. With these values, a study is performed according to the up to date SANTE guide (SANTE/12682/2019 *Guidance document on analytical quality control*), referred to analysis under repeatability conditions. The acceptance criteria to ensure the samples have been stable during the whole proficiency test are the following:

$$\begin{aligned} |(X_{t1} - X_{t2}) / X_{t1}| \cdot 100 &\leq 10\% \\ |(X_{t1} - X_{t3}) / X_{t1}| \cdot 100 &\leq 10\% \end{aligned}$$

Being $|(X_{t1} - X_{tn}) / X_{t1}|$ the difference between the average of the samples analysed before, during and at the end of the proficiency test.

10. EVALUATION REPORT

Once received and statistically evaluated all of the participating laboratories results, TestQual will send a final report that summarizes the participation of each laboratory.

This final report will be received by the laboratories via e-mail in PDF format, but also can be downloaded from the private area of each participant in www.TestQual.com.

If desired, the laboratory may request the report in paper, and it will be sent to its laboratory by mail.

In the event that a participant wishes to appeal against the assessment program performance, a written appellation must be sent by e-mail to jpnavarro@testqual.com explaining the reasons for it.

11. CONTACT

TestQual puts at your disposal any of the following means to contact our team:

Website:	Contact
Email:	jpnavarro@testqual.com
Office phone:	+34 868 94 94 86
Mobile phone:	+34 676 367 555

12. REFERENCES

TestQual Proficiency Testing Schemes are based on the following standards:

UNE-EN ISO/IEC 17043, first edition 2010-02-01. Conformity assessment- General requirements for proficiency testing.

ISO13528:2015, second edition 2015-08-01. Statistical methods for use in proficiency testing by interlaboratory comparison.

THE INTERNATIONAL HARMONIZED PROTOCOL FOR THE PROFICIENCY TESTING OF ANALYTICAL CHEMISTRY LABORATORIES

SANTE/12682/2019, 1st January 2020, Guidance document on analytical quality control and method validation procedures for pesticides residues analysis in food and feed.